



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz column n. Hempel, filling height 200 mm,
sleeve/core NS 14/23, with pushed-on glass
jacket**

**€86.64
plus VAT &
Shipping**

Product Images



Description

Columns according to Hempel

With NS sleeve and NS core, with pushed-on glass jacket.

For filling with packing (e.g. glass beads).

Made of DURAN® tube.

Additional Information

No.	LZ-1056
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3056
EAN	4058072307967
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Form	straight
Height	200 mm
Cooler features	with standard ground joint
Material	Borosilicate glass
Standard ground sleeve	14/23
Standard ground core	14/23

