



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz column head with valve tap, column NS 29/32, template NS 29/32

€979.75
plus VAT &
Shipping

Product Images



Description

Column heads with valve stopcock

With PTFE valve as reflux regulator for easy and accurate reflux adjustment, with ground joint sleeve NS 14/23 for thermometers with an installation length of 55 mm.

Additional Information

| | |
|------------------------|----------------------------|
| No. | LZ-1096 |
| Manufacturer (Brand) | Lenz |
| old neoLab article no. | B-3098 |
| EAN | 4058072308360 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (clear) |
| Material | Borosilicate glass DURAN® |
| Form | square |
| Cooler features | with standard ground joint |
| Material | Borosilicate glass |
| Standard ground sleeve | 29/32 |
| Standard ground core | 29/32 |
| Type cock / valve | PTFE - tap |
| Type cooler | Liebig cooler |

