



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz spherical cooler according to Allihn with
olives, sleeve/core NS 19/26, shell length 160
mm, cooling surface 100 cm²**

**€89.65
plus VAT &
Shipping**

Product Images



Description

Ball cooler according to Allihn with standard ground joint for the condensation of vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: With glass tubing olives. Tubes with an inner diameter of 9–12 mm are required for this.

Additional Information

No.	LZ-1132
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3134
EAN	4058072308728
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 160
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	100 cm ²
Material olives	Borosilicate glass
Standard ground sleeve	19/26
Standard ground core	19/26
Type cooler	Ball cooler
standardized according to DIN standard (DIN)	DIN 12581

