



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




**Lenz ball cooler according to Allihn with GL  
connections, sleeve/core NS 24/29, shell length  
160mm, cooling surface 100 cm<sup>2</sup>**

**€82.35  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Ball cooler according to Allihn with standard ground joint for condensation of vapors in chemical synthesis These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: Water connection: GL 14 thread (with two pre-mounted, unscrewable plastic olives)

## Additional Information

---

No.	LZ-1137
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3139
EAN	4058072308773
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 160
Cooler features	with unscrewable olives with standard grinding
Material	Borosilicate glass
Connection cooling medium	Olive with union nut GL 14
Cooling surface	100 cm <sup>2</sup>
Standard ground sleeve	24/29
Standard ground core	24/29
Type cooler	Ball cooler
Material olives	Polypropylene (PP)

