



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Lenz spherical cooler according to Allihn with  
olives, sleeve/core NS 29/32, shell length 250  
mm, cooling surface 200 cm<sup>2</sup>**

**€72.60  
plus VAT &  
Shipping**

## Product Images

---

## Description

---

Ball cooler according to Allihn with standard ground joint for the condensation of vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: With glass tubing olives. Tubes with an inner diameter of 9–12 mm are required for this.

## Additional Information

---

No.	LZ-1142
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3144
EAN	4058072308827
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 250
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	200 cm <sup>2</sup>
Standard ground sleeve	29/32
Standard ground core	29/32
Type cooler	Ball cooler
Material olives	Borosilicate glass
standardized according to DIN standard (DIN)	DIN 12581

