



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




**Lenz ball cooler according to Allihn with GL connections, sleeve/core NS 24/29, jacket length 250 mm, cooling surface 200 cm<sup>2</sup>**

**€82.84  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Ball cooler according to Allihn with standard ground joint for condensation of vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: Water connection: GL 14 thread (with two pre-mounted, unscrewable plastic olives)

## Additional Information

---

No.	LZ-1145
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3147
EAN	4058072308858
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 250
Cooler features	with unscrewable olives with standard grinding
Material	Borosilicate glass
Connection cooling medium	Olive with union nut GL 14
Cooling surface	200 cm <sup>2</sup>
Material olives	Polypropylene (PP)
Standard ground sleeve	24/29
Standard ground core	24/29
Type cooler	Ball cooler

