



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz spherical cooler according to Allihn with
olives, sleeve/core NS 24/29, shell length 400
mm, cooling surface 300 cm²**

**€82.48
plus VAT &
Shipping**

Product Images



Description

Ball cooler according to Allihn with standard ground joint for the condensation of vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: With glass tubing olives. Tubes with an inner diameter of 9–12 mm are required for this.

Additional Information

No.	LZ-1147
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3149
EAN	4058072308872
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 400
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	300 cm ²
Standard ground sleeve	24/29
Standard ground core	24/29
Type cooler	Ball cooler
Material olives	Borosilicate glass

