



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz ball cooler acc. to Allihn with GL
connections, sleeve/core NS 29/32, shell length
400 mm, cooling surface 300 cm²**

**€95.12
plus VAT &
Shipping**

Product Images

Description

Ball cooler according to Allihn with standard ground joint for condensation of vapors in chemical synthesis These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: Water connection: GL 14 thread (with two pre-mounted, unscrewable plastic olives)

Additional Information

No.	LZ-1150
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3152
EAN	4058072308902
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 400
Cooler features	with unscrewable olives with standard grinding
Material	Borosilicate glass
Connection cooling medium	Olive with union nut GL 14
Cooling surface	300 cm ²
Standard ground sleeve	29/32
Standard ground core	29/32
Type cooler	Ball cooler
Material olives	Polypropylene (PP)

