



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz spherical cooler acc. to Allihn for  
extractors, with olives, core NS 29/32, extractor  
30 ml, cooling surface 100 cm<sup>2</sup>**

**€85.70  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

### Ball cooler according to Allihn with standard ground joint for condensation of vapors in chemical synthesis

These high-quality reflux **condenser** made of borosilicate glass 3.3 are indispensable in the chemical synthesis laboratory nowadays.

#### Design:

Cooler for extractors

With water connection: glass olives or GL 14 thread with plastic olives available

## Additional Information

---

No.	LZ-1153
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3155
EAN	4058072308933
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	30 ml
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	100 cm <sup>2</sup>
Standard ground core	29/32
Type cooler	Ball cooler
Material olives	Borosilicate glass

