



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz ball cooler acc. to Allihn for extractors, with
olives, core NS 45/40, extractor 100 - 250 ml,
cooling surface 200 cm²**

**€99.85
plus VAT &
Shipping**

Product Images



Description

Ball cooler according to Allihn with standard ground joint for condensation of vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: Cooler for extractors. With water connection: glass olives or GL 14 thread available with plastic olives.

Additional Information

No.	LZ-1155
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3157
EAN	4058072308957
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	250 ml
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	200 cm ²
Standard ground core	45/40
Type cooler	Ball cooler
Material olives	Borosilicate glass

