



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




**Lenz ball cooler acc. to Allihn for extractors, with  
GL connections, core NS 29/32, extractor 30 ml,  
cooling surface 100 cm<sup>2</sup>**

**€99.82  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Ball cooler according to Allihn with standard ground joint for condensation of vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: Cooler for extractors. With water connection: glass olives or GL 14 thread available with plastic olives.

## Additional Information

---

No.	LZ-1156
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3158
EAN	4058072308964
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	30 ml
Cooler features	with unscrewable olives with standard grinding
Material	Borosilicate glass
Connection cooling medium	Plastic olive GL 14
Cooling surface	100 cm <sup>2</sup>
Standard ground core	29/32
Type cooler	Ball cooler
Material olives	Polypropylene (PP)

