



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz spherical cooler acc. to Allihn for
extractors, with GL connections, core NS 34/35,
extractor 70 ml, cooling surface 100 cm²**

**€108.04
plus VAT &
Shipping**

Product Images



Description

Ball cooler according to Allihn with standard ground joint for condensation of vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Design: Cooler for extractors. With water connection: glass olives or GL 14 thread available with plastic olives.

Additional Information

No.	LZ-1157
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3159
EAN	4058072308971
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	70 ml
Cooler features	with unscrewable olives with standard grinding
Material	Borosilicate glass
Connection cooling medium	Plastic olive GL 14
Cooling surface	100 cm ²
Material olives	Polypropylene (PP)
Standard ground core	34/35
Type cooler	Ball cooler

