



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz cooler acc. to Liebig with olives,  
sleeve/core NS 24/29, jacket length 250 mm,  
cooling surface 130 cm<sup>2</sup>**

**€50.40  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Standard ground joint reflux condenser for the condensation of gases and vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are indispensable in today's chemical synthesis laboratory. According to DIN 12 756. Liebig coolers: With glass tubing olives. Tubes with an inner diameter of 9–12 mm are used for this purpose.

## Additional Information

---

No.	LZ-1105
Manufacturer (Brand)	Lenz
old neolab article no.	B-3107
EAN	4058072308452
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 250
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	130 cm <sup>2</sup>
Standard ground sleeve	24/29
Standard ground core	24/29
Type cooler	Liebig cooler
Material olives	Borosilicate glass

