



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Lenz cooler acc. to Liebig with GL connections,  
sleeve/core NS 29/32, jacket length 250 mm,  
cooling surface 130 cm<sup>2</sup>**

**€57.05  
plus VAT &  
Shipping**

## Product Images

---

## Description

---

Standard ground joint reflux condenser for condensation of gases and vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are indispensable in today's chemical synthesis laboratory. According to DIN 12 756. Cooler according to Liebig: water connection with thread GL 14, complete with two unscrewable plastic olives.

## Additional Information

---

No.	LZ-1110
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3112
EAN	4058072308506
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 250
Cooler features	with unscrewable olives with standard grinding
Material	Borosilicate glass
Connection cooling medium	Olive with union nut GL 14
Cooling surface	130 cm <sup>2</sup>
Material olives	Polypropylene (PP)
Standard ground sleeve	29/32
Standard ground core	29/32
Type cooler	Liebig cooler
standardized according to DIN standard (DIN)	DIN 12756

