



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz cooler acc. to Liebig with olives,
sleeve/core NS 29/32, jacket length 400 mm,
cooling surface 210 cm²**

**€59.93
plus VAT &
Shipping**

Product Images

Description

Standard ground joint reflux condenser for the condensation of gases and vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are indispensable in today's chemical synthesis laboratory. According to DIN 12 756. Liebig coolers: With glass tubing olives. Tubes with an inner diameter of 9–12 mm are used for this purpose.

Additional Information

No.	LZ-1114
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3116
EAN	4058072308544
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 400
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	210 cm ²
Standard ground sleeve	29/32
Standard ground core	29/32
Type cooler	Liebig cooler
Material olives	Borosilicate glass
standardized according to DIN standard (DIN)	DIN 12576

