



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz cooler acc. to Liebig with GL connections,
sleeve/core NS 29/32, jacket length 400 mm,
cooling surface 210 cm²**

**€64.09
plus VAT &
Shipping**

Product Images

Description

Standard ground joint reflux condenser for condensation of gases and vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are indispensable in today's chemical synthesis laboratory. According to DIN 12 756. Cooler according to Liebig: water connection with thread GL 14, complete with two unscrewable plastic olives.

Additional Information

No.	LZ-1118
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3120
EAN	4058072308582
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 400
Cooler features	with unscrewable olives with standard grinding
Material	Borosilicate glass
Connection cooling medium	Olive with union nut GL 14
Cooling surface	210 cm ²
Material olives	Polypropylene (PP)
Standard ground sleeve	29/32
Standard ground core	29/32
Type cooler	Liebig cooler
standardized according to DIN standard (DIN)	DIN 12756

