



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz snake cooler with olives, sleeve/core NS
14/23, jacket length 160 mm, cooling surface 230
cm²**

**€79.28
plus VAT &
Shipping**



Description

Standard ground joint reflux condenser for condensation of gases and vapors in chemical synthesis

These high-quality reflux condenser made of borosilicate glass 3.3 are indispensable in the chemical synthesis laboratory nowadays.

Snake condenser:

With glass hose olives. Suitable for tubing with an inner diameter of 9-12 mm.

Additional Information

No.	LZ-1159
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3161
EAN	4058072308995
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 160
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	230 cm ²
Standard ground sleeve	14/23
Standard ground core	14/23
Type cooler	Snake cooler
Material olives	Borosilicate glass

