



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Lenz snake cooler with olives, sleeve/core NS  
29/32, jacket length 400 mm, cooling surface  
820 cm<sup>2</sup>**

**€134.02  
plus VAT &  
Shipping**



## Description

---

Standard ground joint reflux condenser for condensation of gases and vapors in chemical synthesis. These high-quality reflux condensers made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Snake cooler: With glass tubing olives. Suitable for tubing with an inner diameter of 9–12 mm.

## Additional Information

---

No.	LZ-1163
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3165
EAN	4058072309022
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 400
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Glass olive
Cooling surface	820 cm <sup>2</sup>
Material olives	Borosilicate glass
Standard ground sleeve	29/32
Standard ground core	29/32
Type cooler	Snake cooler

