



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz serpentine cooler with GL connections,
sleeve/core NS 29/32, shell length 400 mm,
cooling surface 820 cm²**

**€139.48
plus VAT &
Shipping**

Product Images



Description

Standard ground joint reflux condenser for condensation of gases and vapors in chemical synthesis These high-quality reflux condenser made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Coil chiller: Water connection: with GL 14 thread (with two pre-mounted, unscrewable plastic olives)

Additional Information

No.	LZ-1164
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3166
EAN	4058072309039
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 400
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Plastic olive GL 14
Cooling surface	820 cm ²
Standard ground sleeve	29/32
Standard ground core	29/32
Type cooler	Snake cooler
Material olives	Borosilicate glass

