



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




---

**Lenz coil cooler with GL connections,  
sleeve/core NS 29/32, shell length 500 mm,  
cooling surface 1040 cm<sup>2</sup>**

**€159.21  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Standard ground joint reflux condenser for condensation of gases and vapors in chemical synthesis These high-quality reflux condenser made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. Coil chiller: Water connection: with GL 14 thread (with two pre-mounted, unscrewable plastic olives)

## Additional Information

---

No.	LZ-1166
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3168
EAN	4058072309053
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Length	0 500
Cooler features	with standard ground joint
Material	Borosilicate glass
Connection cooling medium	Plastic olive GL 14
Cooling surface	1040 cm <sup>2</sup>
Standard ground sleeve	29/32
Standard ground core	29/32
Type cooler	Snake cooler
Material olives	Borosilicate glass

