



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz cooler with template according to
Twisselmann, core NS 45/40, for extractor 100 -
150 ml**

**€546.88
plus VAT &
Shipping**

Product Images



Description

Cooler with template according to Twisselmann with standard ground joint for the condensation of gases and vapors in chemical synthesis. These high-quality reflux condenser made of borosilicate glass 3.3 are nowadays indispensable in the chemical synthesis laboratory. With 2 pcs. plastic olives GL 14 thread (pre-assembled, can be unscrewed)

Additional Information

No.	LZ-1369
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3372
EAN	4058072311087
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	150 ml
Cooler features	with unscrewable olives with standard grinding
Material	Borosilicate glass
Connection cooling medium	Plastic olive GL 14
Standard ground core	45/40
Type cooler	Twisselmann cooler
Material olives	Polypropylene (PP)

