



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® Schlenkline standard package with Wiggens rotary vane pump, R-8D

€3,684.00
plus VAT &
Shipping

Product Images



Description

The **neoLab Schlenkline standard package** of high-quality individual parts – laced together for your chemical synthesis laboratory. With powerful rotary vane pumps, a high vacuum of 0.002 mbar can be achieved. This ensures permanent water- and air-free synthesis. The Schlenkline – handmade by experienced glassblowers in Germany – is made of durable borosilicate glass, 5 patent 2-way glass faucets and olives for tubing with an inner diameter of 9 mm.

Scope of delivery:

- Wiggins rotary vane pump F-1003
- Schlenkline
- Cold trap 100 ml
- suitable Dewar vessel
- 5 m vacuum hose, red
- 5 m silicone hose for nitrogen/argon, transparent
- 5 stainless steel hose clamps 12-22 mm
- 5 stainless steel hose clamps 16-27 mm
- 2 gas wash bottles
- Three-way valve made of PP with PTFE plug

All parts also available separately at [neolab.de](https://www.neolab.de)

Additional Information

No.	E-0001
Manufacturer (Brand)	neoLab
EAN	4058072160937
Transport temperature	Room temperature
Material	Borosilicate glass

