



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Combitips advanced Adapter 50 ml, light gray

**€11.80**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Adapter for inserting the Combitips advanced® into the Multipette®. Suitable for Multipette® plus, Multipette® M4, Multipette stream® or Xstream®. For Combitips advanced with a volume of 50 ml.

## Additional Information

---

No.	VB-0265
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-6410
EAN	4058072405458
Transport temperature	Room temperature

