



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex PP Adapter for Calibrex™ 520 Dispenser, thread Ø25 mm

**€8.84**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The Calibrex® dispensers are supplied with 3 thread adapters each.  
For other bottle diameters, select an appropriate adapter or adapter combinations to obtain the desired diameter.

The 520 model has a 32 mm base thread.

## Additional Information

---

No.	SR-0185
Manufacturer (Brand)	Socorex
old neolab article no.	7-7457
EAN	4058072367039
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
suitable for	Bottles, round
DM outside	25 mm
Color	Gray
Material	Polypropylene (PP)
Thread	GL 25
suitable for	Fixing/ Connecting

