



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex PP adapter for Calibrex™ 520 dispenser, thread Ø30 mm

**€17.68**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The Calibrex® dispensers are supplied with 3 thread adapters each.  
For other bottle diameters, select an appropriate adapter or adapter combinations to obtain the desired diameter.

The 520 model has a 32 mm base thread.

## Additional Information

---

|                        |                    |
|------------------------|--------------------|
| No.                    | SR-0187            |
| Manufacturer (Brand)   | Socorex            |
| old neolab article no. | 7-7459             |
| EAN                    | 4058072367053      |
| Transport temperature  | Room temperature   |
| Color                  | Black              |
| Material               | Polypropylene (PP) |
| suitable for           | Bottles, round     |
| DM outside             | 30 mm              |
| Color                  | Gray               |
| Material               | Polypropylene (PP) |
| Thread                 | Special thread     |
| suitable for           | Fixing/ Connecting |

