



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex PP Adapter for Calibrex™ 520 Dispenser, Thread Ø36 mm

€26.94
plus VAT &
Shipping

Product Images



Description

The Calibrex® dispensers are supplied with 3 thread adapters each.
For other bottle diameters, select an appropriate adapter or adapter combinations to obtain the desired diameter.

The 520 model has a 32 mm base thread.

Additional Information

No.	SR-0189
Manufacturer (Brand)	Socorex
old neolab article no.	7-7461
EAN	4058072367077
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
suitable for	Bottles, round
DM outside	36 mm
Color	Gray
Material	Polypropylene (PP)
Thread	GL 36
suitable for	Fixing/ Connecting

