



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex PTFE adapter for Calibrex™ 520 dispenser, thread Ø38 mm

**€79.46**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The Calibrex® dispensers are supplied with 3 thread adapters each.  
For other bottle diameters, select an appropriate adapter or adapter combinations to obtain the desired diameter.

The 520 model has a 32 mm base thread.

## Additional Information

---

No.	SR-0194
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7466
EAN	4058072367121
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
suitable for	Bottles, round
DM outside	38 mm
Color	Gray
Material	Polypropylene (PP)
Thread	GL 38
suitable for	Fixing/ Connecting

