



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex PP Adapter Ø 32-28 mm for Calibrex™ organo 525 and solutae 530 Dispensers

€17.58
plus VAT &
Shipping

Product Images



Description

The Calibrex® dispensers are supplied with 3 thread adapters each.
For other bottle diameters, select an appropriate adapter or adapter combinations to obtain the desired diameter.
Models 525 and 530 have a 45 mm base thread.

Additional Information

No.	SR-0233
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7610
EAN	4058072367510
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
suitable for	Bottles. round
DM outside	32 mm
DM inside	28 mm
Color	Gray
Material	Polypropylene (PP)
Thread	GL 45
suitable for	Fixing/ Connecting

