



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN C-mount camera adapter 0.60x (adjustable focus)

**€175.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- C-mount camera adapter 0.60× (adjustable focus)

## Additional Information

---

No.	KP-4781
Manufacturer (Brand)	KERN
EAN	4058072143718
Transport temperature	Room temperature
Type Connection Input	C-Mount
Manufacturer	CORE

