



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine silicone adapter for D-6018 and D-6020

€8.54
plus VAT &
Shipping

Product Images



Description

Silicone adapter suitable for pipettes with a volume of 0.1 to 100 ml.
Autoclavable.

Additional Information

No.	D-6026
Manufacturer (Brand)	neoLabLine
EAN	4058072118921
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	16 mm
Color	transparent (milky)
Material	Silicone
Height	53 mm
DM Bore	6 mm
for model pipette controller	neoLab Line pipette controller. yellow neoLab Line pipette controller. pink

