



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



sunlab® replacement silicone adapter for D-8701

€10.30
plus VAT &
Shipping

Product Images



Description

Replacement silicone adapter for Sunlab pipette controller SU1700 (D-8700) Suitable for 0.1-100 ml pipettes; autoclavable at 121°C PU: 1 piece

Additional Information

No.	D-8715
Manufacturer (Brand)	sunlab
EAN	4058072088385
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Silicone
for model pipette controller	Sunlab® digital pipette controller SU1700

