




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Paste pipette adapter for Acura® 835.02 and Calibra® 832.03

€29.74
plus VAT &
Shipping

Product Images



Description

Adapter for Pasteur pipettes made of PVDF

With 2 internal O-rings for a tight fit of the glass Pasteur pipettes (AD 6.5–7.2 mm).

Additional Information

No.	SR-0408
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7857
EAN	4058072369217
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Silicone
for model pipette controller	Socorex Acura® 832.02 Calibra® 832.02

