




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bohlender b.safe septum adapter UNF 1/4", 1x ND9

€32.38
plus VAT &
Shipping

Product Images



Description

b.safe Septum Adapter Threaded adapter for connecting a septum cap to the capillary connection of b.safe Caps.
Dimensions: Suitable for capillary connector / Capillary connector (UNF 1/4) Transition to septum cap / Septum cap (ND9) Material PTFE "

Additional Information

No.	ER-2601
Manufacturer (Brand)	Bohlender
EAN	4058072526948
Transport temperature	Room temperature
Net weight	0.00390 kg
Country of origin	Germany
Customs tariff number	39269097

