



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® One-hand shaking attachment,
rubber. Support, for vessels up to 100 ml, for
3-5096**

**€53.71
plus VAT &
Shipping**

Product Images



Description

Order No. 3-5099 One-hand shaking attachment Fig. 4

This shaking attachment with rubberized support is suitable for simultaneous mixing of several test tubes or flasks up to 100 ml. The utensils must be held by hand during the shaking process.

Dimensions: 137 x 132 x 44 mm (L x W x D)

Recommended max. shaking frequency: 1000 rpm

Additional Information

No.	ZT-0572
Manufacturer (Brand)	neoLab
old neoLab article no.	3-5099
EAN	4058072629274
Transport temperature	Room temperature
Figure	4
Attachment is suitable for type vessel	Test tubes flasks
for model shaker/ rotator	neoLab universal shaker (3-5095)

