



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Shaker attachment for 4 flasks up to 500 ml, for shaker 3-5096

**€131.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### **Order No. 3-5097 Bottle shaker attachment Fig. 2**

Suitable for up to 4 **round-bottomed** or Erlenmeyer flasks up to 500 ml capacity each. Due to the low center of gravity, particularly smooth running is ensured. **Clamp span:** 60 mm

Recommended max. shaking frequency: 500 rpm

## Additional Information

---

No.	ZT-0574
Manufacturer (Brand)	neoLab
old neoLab article no.	3-5097
EAN	4058072629298
Transport temperature	Room temperature
Figure	2
Attachment is suitable for type vessel	Piston
for model shaker/ rotator	neoLab universal shaker (3-5095)

