



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine turntable for 20 tubes 8-12 mm Ø, for rotator

€132.60
plus VAT &
Shipping

Product Images



Description

The angle of the **turntable** can be adjusted between 27° and 90° and can be replaced in just a few steps. A suitable adapter (2-1183) for using **two turntables** can be found in the accessories. Vessels are held firmly by the vinyl-covered spring steel clamps.

Further turntables on request.

Additional Information

No.	2-1179
Manufacturer (Brand)	neoLabLine
EAN	4058072055752
Transport temperature	Room temperature
for vessels with volume up to	10 ml
Spare part	1-2131 1-2132
Material turntable	Aluminum
Material clamp	Stee powder coated
for vessels with DM up to	12.0 mm
for model shaker/ rotator	neoLabLine Rotator 10 rpm fixed (2-1175) neoLabLine Rotator 5 to 50 rpm (2-1184) neoLabLine Rotator 5 to 50 rpm digital (2-1185)
for type (micro) vessel	Reaction vessels

