



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## sunlab® plate with clamps 60 x for 1.5 ml reaction tubes

**€173.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Turntable for Sunlab Rotator, digital

The plate is fixed on the rotator.

A suitable adapter (D-8504) for using **two turntables** can be found in the accessories.

It is equipped with clamps in which reaction tubes and centrifuge tubes of different volumes can be clamped.

**Order No. D-8501:** 60 x 1.5 ml reaction tubes

## Additional Information

---

No.	D-8501
Manufacturer (Brand)	sunlab
EAN	4058072024376
Transport temperature	Room temperature
for vessels with volume up to	1.5 ml
Spare part	D-8505
Material turntable	Aluminum
Material clamp	Steel
for vessels with DM up to	11.0 mm
for model shaker/ rotator	Sunlab® Rotator. digital (D-8500)
for type (micro) vessel	Reaction vessels

