



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



sunlab® plate with clamps 8 x for 50 ml tubes

€131.90
plus VAT &
Shipping

Product Images



Description

Turntable for Sunlab Rotator, digital

The plate is fixed on the rotator.

A suitable adapter (D-8504) for using **two turntables** can be found in the accessories.

It is equipped with clamps in which reaction tubes and centrifuge tubes of different volumes can be clamped.

Order No. D-8503: 8 x 50 ml centrifuge tubes

Additional Information

No.	D-8503
Manufacturer (Brand)	sunlab
EAN	4058072024390
Transport temperature	Room temperature
for vessels with volume up to	50 ml
Spare part	D-8507
Material turntable	Aluminum
Material clamp	Steel
for vessels with DM up to	28.0 mm
for model shaker/ rotator	Sunlab® Rotator. digital (D-8500)
for type (micro) vessel	Centrifuge tube

