



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



sunlab® Replacement Adapter from 2 ml to 0.2 ml PCR

€23.10
plus VAT &
Shipping

Product Images



Description

Replacement inserts for Rotor Sunlab Mini Centrifuge

Reducing insert for centrifugation of small-volume tubes.

Order No. D-8551: 2 ml to 0.2 ml (PCR tubes)

Additional Information

No.	D-8551
Manufacturer (Brand)	sunlab
EAN	4058072024437
Transport temperature	Room temperature
Adapter is required for	PCR tubes. 8 x 0.2 ml (D-8551)
for model centrifuge	Sunlab Mini Centrifuge

