



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



sunlab® flask attachment for 500 ml

€140.60
plus VAT &
Shipping

Product Images



Description

Attachments for Sunlab magnetic stirrers

Made of aluminum, suitable for round flasks.

Capacity: 1 flask each

Please note:

When the flask attachments are in use, please do not leave the magnetic stirrer unattended !

The glass flasks must be removed from the flask attachment during the cooling process !

Due to the different cooling rates of glass and aluminum, the glass vessel may otherwise break.

Additional Information

No.	D-8152
Manufacturer (Brand)	sunlab
EAN	4058072018030
Transport temperature	Room temperature
Material	Aluminum
Attachment is suitable for	Round bottom flask 500 ml (D-8152)
for vessels with volume up to	500 ml

