



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® polycarbonate mixing attachment 4 l, for rotor blender

Special Price
€250.00 was
€375.00
plus VAT &
Shipping

Product Images



Description

Mixer-homogenizer for chopping, mixing, stirring, emulsifying and homogenizing.

Mixing attachments with special blade heads and specially shaped blades made of stainless special blade steel. Low center of gravity and maximum stability due to extra heavy motor housing made of chrome-plated die-cast zinc. Wear-resistant and vibration-free direct drive through collector motor.

The cups have an anatomically shaped insulated handle for one-handed decanting without dripping. The base pieces of the cups are made of metal, easily removable by a small twist for hygienic cleaning.

- Continuous power 800 W, peak power up to 1000 W
- infinitely variable speed from 600 to 14,000 rpm
- Soft start
- Motor protection switch with electronic overload and overheating monitoring

Safety control:

The motor only runs when the cover with retaining bar is properly fitted.

CE-TÜV-GS

Connection: 230 Volt

Weight: 6 kg

Height: approx. 50 cm

Additional Information

No.	3-2067
Manufacturer (Brand)	neoLab
EAN	4058072018764
Transport temperature	Room temperature
Color	transparent (clear)
Type closure/ lid	Plug closure (inside)
Color graduation	transparent (milky)
Division	200 ml

