



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® 2 liter mixing attachment poycarbonate

€312.00
plus VAT &
Shipping

Product Images



Description

Mixing attachments for rotor high-performance laboratory mixer, programmable 7 different mixing attachments: 125 ml, 250 ml, 500 ml, 1000 ml, 2000 ml The mixing attachments 125 to 1000 ml are closed at the top and are made of glass. The 2000 ml mixing tips are made of polycarbonate or stainless steel and have baffles for better mixing.

Additional Information

No.	3-2045
Manufacturer (Brand)	neoLab
EAN	4058072018641
Transport temperature	Room temperature
Color	transparent (milky)
Type closure/ lid	Overlay closure
Color graduation	transparent (milky)
Division	200 ml

