



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



sunlab® neoLab Additional Rack f. SU1801, 1.5 ml, vertical/horizontal shaking motion

€103.00
plus VAT &
Shipping

Product Images



Description

The various attachments are suitable for reaction tubes and centrifuge tubes from 1.5 to 50 ml. Depending on the attachment, the tubes can be locked vertically, horizontally or variably in the rotor rack.

Additional Information

No.	D-8802
Manufacturer (Brand)	sunlab
EAN	4058072085537
Transport temperature	Room temperature
Material	Stainless Steel Acrylonitrile Butadiene Styrene (ABS)
Attachment is suitable for	Reaction tubes. 36 x 1.5/2.0 ml (D-8802)
Attachment is suitable for type vessel	Reaction vessels
Layer tube	vertical
for vessels with DM up to	12.5 mm

