



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



sunlab® neolab additional rack for SU1801, 5/ 15 ml, horizontal shaking movement

€103.00
plus VAT &
Shipping

Product Images



Description

The various attachments are suitable for reaction tubes and centrifuge tubes from 1.5 to 50 ml. Depending on the attachment, the tubes can be locked vertically, horizontally or variably in the rotor rack.

Additional Information

No.	D-8805
Manufacturer (Brand)	sunlab
EAN	4058072085568
Transport temperature	Room temperature
Material	Stainless Steel Acrylonitrile Butadiene Styrene (ABS)
Attachment is suitable for	Reaction vessels and centrifuge tubes. 12 x 5.0/15 ml (D-8805)
Attachment is suitable for type vessel	Reaction vessels Centrifuge tubes
Layer tube	horizontal
for vessels with DM up to	17.0 mm

