



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Aluminum drying bath, PTFE jacket

€438.30
plus VAT &
Shipping

Product Images



Description

The **neoLabLine drying bath made of aluminum** and a PTFE jacket offers the pleasant alternative to an oil bath. Clean, simple and uncomplicated, yet just as efficient as an oil bath.

Simply fill the drying bath with an appropriate amount of heating beads, and then place the vessel with your sample in it. The desired temperature is quickly reached due to the good thermal conductivity of the special metal alloy.

- No more cleaning, changing and refurbishing necessary
- Cost-saving
- Environmentally friendly

Additional Information

No.	D-6111
Manufacturer (Brand)	neoLab
EAN	4066292501010
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Aluminum Polytetrafluoroethylene (PTFE. Teflon®)
Attachment is suitable for	different vessels
for vessels with volume up to	2000 ml

