



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Replacement pouring ring, 16 pcs./pack

€11.04
plus VAT &
Shipping

Product Images



Description

PP pouring rings for ^{DURAN®} YOUTILITY bottles

pouring ring GL 45

Color: cyan

PU: 16 pieces

Additional Information

No.	SD-2755
Manufacturer (Brand)	DURAN
old neoLab article no.	2-3048
EAN	4058072363604
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Polypropylene (PP)
Height	3 mm
DM inside	45 mm
TBST MAX pouring ring	140 C°
TBST MIN pouring ring	5 C°

