



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® Original pouring ring, for laboratory bottles, PP, yellow

**€6.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

DURAN® Original GL 45 pouring ring, for laboratory bottles, PP, yellow Available in blue, green, yellow and gray with matching color pouring rings. Keeping media apart is simplified, interchanging screw caps and carry-over of substances is virtually eliminated. Special features and benefits: With integrated lip seal and pouring ring for drip-free removal and clean, safe working With pouring ring for drip-free removal Clear labeling prevents interchanging of closures and thus carryover of substances Compatible with all GL 45 threads and thus versatile Application: Autoclaving of media and reagents, short-term storage of samples and media, color labeling and identification

## Additional Information

---

No.	SD-0015
Manufacturer (Brand)	DURAN
VGKL number	112989917
EAN	4058072187019
Manufacturer EAN	4032051246867
Gross weight	0.001
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.00100 kg
Color	yellow
Material	Polypropylene (PP)
Country of origin	Germany
sterile	No
Region of origin	Rhineland-Palatinate
Customs tariff number	39269097

