



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® silicone identification rings, 16 pcs./pack

€29.28
plus VAT &
Shipping

Product Images



Description

DURAN® Bottle Tag GL 45

The colored elastic identification rings made of silicone are used for personal labeling and clear identification. They can be easily and safely applied to the neck of the bottle, even if the bottle is already fitted with the screw cap.

2 pcs. per color blue, purple, green, red, orange, yellow, black, white for GL 45 bottles.

Additional Information

No.	SD-2766
Manufacturer (Brand)	DURAN
old neoLab article no.	2-3042
EAN	4058072363703
Transport temperature	Room temperature
Color	blue purple green red orange yellow black white
Material	Silicone
Number of parts in set	16

