




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® button cell V13GA (LR44), alkaline, 1.5 V, 125 mAh

€4.88
plus VAT &
Shipping

Product Images



Description

High performance button cells for a wide variety of electronic devices.

Button cells supplied by us can be returned to us for disposal.

Additional Information

| | |
|---------------------------------------|---|
| No. | 4-1080 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072020729 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| sterile | No |
| Batteries and accumulators properties | not rechargeable |
| Color | Silver |
| Size | Button cell |
| Electrochem. System | primary alkali-manganese (Zn/MnO ₂) |
| Voltage | 1.5 V |
| Type battery | Battery (coin cell) |
| Type | V13GA |

