



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® button cell 4001 (LR1), alkaline, 1.5 V, 800 mAh

€5.83
plus VAT &
Shipping

Product Images



Description

Brand alkaline batteries, mercury and cadmium free.

Very high energy content and good voltage level during discharge. Universally applicable for all devices. These batteries are not rechargeable. Batteries supplied by us can be returned to us for disposal.

Additional Information

No.	4-2081
Manufacturer (Brand)	neoLab
EAN	4058072022730
Transport temperature	Room temperature
autoclavable at 121°C	No
sterile	No
Batteries and accumulators properties	not rechargeable
Color	Silver
Size	Button cell
Electrochem. System	primary alkali-manganese (Zn/MnO ₂)
Voltage	1.5 V
Type battery	Battery (coin cell)

